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DVR 3pro  II 2 G Ex ib IIC T4 Gb

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<b>A</b>		<b>L</b>	
Accessories	43	Limitation of use	18
Action step	8	<b>M</b>	
approved battery types	43	Mandatory sign	7
Approved block-type battery	43	Measuring cycle	26, 29
ATEX specification	15	<b>O</b>	
Atmospheric pressure	33	Operating elements	23
<b>B</b>		Operating elements DVR 3pro	23
Block-type battery (battery list)	43	Ordering information	43
<b>C</b>		<b>P</b>	
captive screw	16	Personnel (staff)	12
Cleaning	31	Presentation convention operating steps	8
Clean sensor	31	Prohibition symbol	7
Clean surface	31	Proper disposal	13
Connection options	20	<b>R</b>	
Contact	5	Rating plate	42
Copyright ©	6	Rating plate data	42
Correct battery handling	13	Rear view	16
<b>D</b>		Replace battery	19
Delivery state automatic switch-off	28	Resolving problems	38
Device data	42	Return (reshipment)	44
Device name	15	<b>S</b>	
Display battery discharge state	30	Safety	5
Display elements DVR 3pro	25	Safety Information for Vacuum Equipment	10
Display icons DVR 3pro	26	Safety instructions	10
<b>E</b>		Safety precautions	12
EC Declaration of Conformity	46	Sales offices	43
Entsorgung	13	Sensor	40
Error display	38	Service range	44
Examples of use:	17	Servicing handling	44
<b>F</b>		Side view	15
Fault – Cause – Remedy	39	Signals	8
Foreseeable misuse	11	Source of supply	43
Front and rear side	15	Spare parts	43
<b>G</b>		Symbols	7
Gauge description	15	<b>T</b>	
Goods arrival	14	Technical data	40, 41
<b>H</b>		Technical information	40
Handling instruction	8	Technical support	38
Health and Safety Clearance	44	<b>U</b>	
<b>I</b>		User information	5
Icons	7	<b>V</b>	
Improper use	10, 11	Vacuum connection	20
Included materials	14	Vacuum sensor	40
Insert battery	19	<b>W</b>	
Installation	18	Warning levels	6
Installation conditions	18	Warning symbol	7
Intended use	10	Wetted materials	41
<b>J</b>			
Jump-back times	24		
<b>K</b>			
Key combination	24		
Key function	23		

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VACUUBRAND GMBH + CO KG has made reasonable efforts to ensure that hazardous materials and substances may not be used in its products.

In order to determine the concentration of hazardous substances in all homogeneous materials of the subassemblies, a “Product Conformity Assessment” (PCA) procedure was performed. As defined in GB/T 26572 the “Maximum Concentration Value” limits (MCV) apply to these restricted substances:

- Lead (Pb): 0.1%
- Mercury (Hg): 0.1%
- Cadmium (Cd): 0.01%
- Hexavalent chromium (Cr(+VI)): 0.1%
- Polybrominated biphenyls (PBB): 0.1%
- Polybrominated diphenyl ether (PBDE): 0.1%

### Environmentally Friendly Use Period (EFUP)

EFUP defines the period in years during which the hazardous substances contained in electrical and electronic products will not leak or mutate under normal operating conditions. During normal use by the user such electrical and electronic products will not result in serious environmental pollution, cause serious bodily injury or damage to the user's assets.

The Environmentally Friendly Use Period for VACUUBRAND products is 40 years.

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This table is created according to SJ/T 11364-2014.

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**Note:** Table applies to all products. Some of the components or parts listed above may not be part of the enclosed product.

- O: 表示该有毒有害物质在该部件所有均质材料中的含量均在GB/T 26572规定的限量要求以下。  
 O: Indicates that the above mentioned hazardous substance contained in all homogeneous materials of the part is below the required limit as defined in GB/T 26572.
- X: 表示该有毒有害物质至少在该部件某一均质材料中的含量超出GB/T 26572规定的限量要求。  
 X: Indicates that the above mentioned hazardous substance contained in at least one of the homogeneous materials of this part is above the required limit as defined in GB/T 26572.

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Apart from the disclosures in the above table, the subassemblies are not intentionally manufactured or formulated with lead (Pb), mercury (Hg), cadmium (Cd), hexavalent chromium (Cr+VI), polybrominated biphenyls (PBB), and polybrominated diphenyl ethers (PBDE).

Products manufactured by VACUUBRAND may enter into further devices (e.g., rotary evaporator) or can be used together with other appliances (e.g., usage as booster pumps).

With these products and appliances in particular, please note the EFUP labeled on these products. VACUUBRAND will not take responsibility for the EFUP of those products and appliances.

Place, date: Wertheim, 06/04/2020

(Dr. F. Gitmans)  
 Managing Director

i.A. (Dr. A. Wollschläger)  
 Regulatory Affairs Manager

**VACUUBRAND GMBH + CO KG**  
 Alfred-Zippe-Str. 4  
 97877 Wertheim  
 Germany

Tel.: +49 9342 808-0  
 Fax: +49 9342 808-5555  
 E-Mail: [info@vacuubrand.com](mailto:info@vacuubrand.com)  
 Web: [www.vacuubrand.com](http://www.vacuubrand.com)

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**VACUUBRAND GMBH + CO KG**  
**Alfred-Zippe-Str.4**  
**97877 Wertheim**  
**德国**

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+49 9342 808-5555

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[info@vacuubrand.com](mailto:info@vacuubrand.com)

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[www.vacuubrand.com](http://www.vacuubrand.com)