



Technology for Vacuum Systems


Have you heard...?

The new remote control
VACUU · CONTROL™



...the power's in your hands! Remote control via VACUU · CONTROL™ !

Available as a modular accessory for all CVC 3000 and DCP 3000
from software version 2.0 onwards (2009).





VACUU · CONTROL™ enables web based remote monitoring and control of vacuum pumping units via computers or mobile devices such as Smartphones and Tablets.

Now it is easy to integrate your vacuum system into your computer network via the [LAN adapter](#) or wireless [WLAN adapter](#) without the need for any additional installation.



REMOTE CONTROL

The graphical user interface allows clear monitoring, documentation and remote control of your vacuum systems.



DOCUMENTATION

The integrated data logger function permanently records the pressure curves of your applications. This simplifies the documentation of your processes.



MONITORING

You can monitor the vacuum processes in your laboratory from the comfort of your office via smartphone, tablet, laptop and PC. So you can take corrective action via the remote software at any time.

ORDERING INFORMATION

VACUU · CONTROL™ WLAN 683110

Remote control adapter WLAN for integration into WLAN PC networks, with web based, graphical user interface, complete with cable for connection to RS 232C interface for CVC 3000 or DCP 3000.

VACUU · CONTROL™ LAN 683120

Remote control adapter LAN with RJ45 socket for integration into LAN PC network, with web based, graphical user interface, complete with cable for connection to RS 232C interface for CVC 3000 or DCP 3000.

VACUUBRAND GMBH + CO KG

Alfred-Zippe-Straße 4
97877 Wertheim · Germany

T +49 9342 808-5550

F +49 9342 808-5555

info@vacuubrand.com

www.vacuubrand.com



Technology for Vacuum Systems