


LOCAL AREA VACUUM NETWORK

VACUU·LAN® - VACUUM FOR MANY USERS



- For new and existing laboratories
- One laboratory, one vacuum pump, several applications and workstations
- Individual configuration with chemical resistant PTFE tubing
- Space-saving installation with small investment
- Demand controlled energy saving vacuum
- Independent work without reciprocal influence - via integrated non-return valves
- Flexibility - easy upgrade and quick replacement of components
- Environmentally friendly - solvent recovery near 100%
- Higher safety for the laboratory - additional components (e.g., emission condenser Peltronic®, flask with liquid level sensor, solvent drain)
- Minimal energy consumption and maximum maintenance intervals by pump on demand control with VARIO® technology

Ask for  
VACUU·LAN® information material