



Changelog vacuum gauge VACUUBRAND DCP 3000

Note: Minor bugfixes not all explicitly listed.

Date	Firmware	Changes
01.10.2015	V2.32	Minor bugfixes
13.05.2015	V2.31	V2.30 issue regarding VSP 3000 adjustment solved
27.02.2015	V2.30	Compatibility with remote control VACUU-CONTROL®
11.02.2013	V2.27	Minor bugfixes
14.11.2013	V2.26	Minor bugfixes
16.04.2012	V2.25	Minor bugfixes
14.03.2012	V2.24	Flexible RS232 input format for pressure settings
30.11.2010	V2.23	Minor bugfixes
01.07.2010	V2.22	Up to 4 MPT 100 gauges heads can be connected
05.05.2010	V2.21	RS 232: Pressure value output with leading zeros RS 232: Output with +- sign for reference measurements
11.03.2010	V2.20	Analog output for every gauge head possible Switching point functionality included MPT 100 output via analog I/O module (0.75 V/decade)
11.09.2009	V2.10	New language: Finnish
22.07.2009	V2.03	Minor bugfixes for MPT 100 gauge head
08.06.2009	V2.02	Manufacturing optimized
20.05.2009	V2.01	Adjustment of MPT 100 gauge head possible
12.05.2009	V2.00	New hardware. Note: Software >2.xy not compatible with older devices 13 new languages MPT 100 gauge head
19.11.2008	V1.31	Factory settings optimized
31.07.2008	V1.30	Data logger Differential pressure indication Graphical pressure over time display
24.04.2008	V1.20	Analog I/O module also for fine vacuum gauge head (logarith. output)
17.12.2007	V1.10	Language selection on first switch on Improvement of VACUU-BUS configuration
22.11.2007	V1.01	Improvement of switching time between gauge heads
08.05.2007	V1.00	First release