

Connecting parts for rotary evaporators and VACUUBRAND vacuum pumps



**MZ 2C NT + AK + EK**

- 7 mbar
- 2.0 m<sup>3</sup>/h

without vacuum control for solvents with boiling point < 100 °C



**MZ 2C NT + AK + EK**

- 2 mbar
- 1.5 m<sup>3</sup>/h

without vacuum control for solvents with boiling point > 100 °C



**MZ 2C NT + AK + EK**

- 2 mbar
- 2.0 m<sup>3</sup>/h

automatic VARIO®-control for solvents with > 100 °C

**Manufacturer**

**Büchi**

- R-300 Serie
- R-200 Serie
  - R-200 / R-205
  - R-210 / R-215
- R-II Serie
- R-100 Serie
  - R-100
  - R-114 / R-124 / RE-121

Vacuum tubing DN 6  
Item-no. 20686000  
Hose nozzle DN 6/10mm  
Item-no. 20636635

Vacuum tubing DN 6  
Item-no. 20686000

**Heidolph**

- Hei-VAP Serie
  - Value / Value digital
  - Advantage / Precision
- Laborata 4000 Serie
- 4000 / 4001 efficient
  - 4010 / 4011 digital
  - 4002 / 4003 control

Vacuum tubing DN 8  
Item-no. 20686001

**IKA**

- RV 8
- RV 10 digital / control

Vacuum tubing DN 8  
Item-no. 20686001

**KNF**

- RC 600
- RC 900

Vacuum tubing DN 8  
Item-no. 20686001